



Bringing you value
and support

BAUSCH+LOMB

Boston[®]
Materials

Bausch + Lomb, Boston® Materials

Adding Value.

Expand the custom contact lens business through product, services, technology

- GP materials
- Soft lens materials
- Licensing of designs
- Educational resources
- Regulatory knowledge
- Industry support
- Strength of B+L and Boston brand to consumers and eye care professionals

Material matters!



Global Leader for Decades

Boston® GP materials

- US indication of ocular surface disease and dry eye
- Customized buttons
- Supported by guarantee
- Hologram stickers available


Soft lens materials

Hydra PEG coating

Educational Resources

Educational guides

- Scleral Lens Guide
- Keratoconus Guide
- Presbyopia Guide
- Multiple Languages
- Online and printed formats



Online Tools:

- Fitting videos
- Webinars



Conferences and Workshops

- Expert speakers



Timeline of Boston® Materials

- 1983: Boston™ II silicone acrylate material is introduced.
- 1983: Boston™ solutions are introduced in the United States.
- 1984: Boston™ IV silicone acrylate material is introduced.
- 1986: Polymer Technology Corporation introduces its first fluoro silicone acrylate material, Boston™ Equalens™.
- 1991: Boston™ Equalens™ II fluoro silicone acrylate lens material is introduced.
- 1994: Boston 7™ AERCOR™ low silicon/high stability fluoro silicone acrylate lens material is introduced.
- 1995: Boston ES™ with AERCOR™ architecture is introduced.
- 1997: Boston XO™ with high Dk/high stability is introduced.
- 1998: Boston EO™ with AERCOR™ architecture is introduced.
- 2005: Bausch & Lomb Vision Shaping Treatment VST™ is introduced for overnight orthokeratology (U.S. market only).
- 2006: Plasma treatment is FDA approved for all Boston materials.
- 2007: Boston XO₂™, a hyper-Dk material, is introduced.
- 2008: Guide to Keratoconus and Guide to Scleral Lenses v 1.0 introduced
- 2013: Scleral Fitting Video series launched
- 2014: Guide to Scleral Lenses v2.0 introduced
- 2017: Boston XO and XO₂™ US indication for ocular surface disease and dry eye
- 2018: HydraPeg Coating Licensing Agreement on Boston ES, EO, XO and XO₂

Boston® Materials

Introduction to Gas Permeable Materials

Bausch + Lomb Boston's commitment to the field of GP contact lens research:


- improve the comfort,
- physiology,
- and safety of GP lens wear.

The following section contains a basic explanation of corneal physiology and a description of various physical characteristics of gas permeable contact lens materials.

This information permits the fitter to select the appropriate gas permeable lens material for each wearer.

Boston® Materials Corneal Physiology

The cornea has 5 layers. They contribute to corneal strength and function.



One of the most important of these corneal layers is the endothelium.

Functions of the endothelium:

- Permeability barrier, allowing the diffusion of nutrients to the cornea.
- Pump mechanism regulates water to maintain the cornea in a partially hydrated state.

Endothelial cells have limited (if any) capacity to regenerate.

Changes in endothelial morphology can occur for example from:

- trauma,
- inflammation,
- long-term oxygen deprivation (hypoxia) as seen in wearers of PMMA.

It is for these reasons that GP lens permeability plays an important role in corneal physiology.

Boston® Materials

Physical Properties of GP Materials

Oxygen

Dk: This term is used to denote the oxygen permeability of rigid and soft contact lens materials.

“**D**” = inherent ability of the material to allow oxygen to diffuse through;

“**k**” = degree to which oxygen is solubilized within or on the material, since water plays some role in absorbing and assisting in the transport of oxygen.

Several methods to measure permeability.

- **Gas-to-gas method** uses a wafer of lens material or a lens affixed to a graduated capillary tube. This allows the volume of oxygen to be measured as it passes through the test material.
- **ISO/Fatt method** uses a wafer or lens affixed to the end of an oxygen probe. The probe is then immersed in a liquid medium. The amount of oxygen that is able to flow through the lens material is then measured. ISO/Fatt is the method used to determine the Dk of Boston materials.

Boston® Materials

Physical Properties of GP Materials


Oxygen

Dk/t: Refers to the transmissibility of a material when it is made into either a plus or minus lens; “t” represents the thickness of a given lens. The significance of this measurement is that the amount of oxygen transmitted can vary depending on the thickness.

EOP: This is perhaps a more meaningful and clinically important value, since it represents the actual amount of oxygen that passes through the lens and reaches the cornea. This measurement takes into account the total lens (material and design).

EOP - Equivalent oxygen percentage

(ISO/Fatt method expressed in cgs units. EOP and Dk/t measurements are approximate. Data on file.)



Boston® Materials

Physical Properties of GP Materials

Wetting angle

Wetting angle is often used as a predictor of the on-eye wetting characteristic of a GP contact lens material.

In theory, a low contact angle equates to good lens surface wetting.

The standard test methods for measuring the contact angle of contact lens materials are defined in ANSI Z80.20-2004 as the sessile drop method and the captive bubble method.

| | PMMA | Boston® II | BostonES® | Boston® IV | Boston® Equalens® | Boston EO® | Boston® Equalens® II | Boston XO® | Boston XO ₂ ® |
|-----------------|------|------------|-----------|------------|----------------------|------------|-------------------------|------------|-----------------------------|
| Captive Bubble* | 60° | 20° | 52° | 17° | 30° | 49° | 30° | 49° | 38° |

* Data on file |

Boston® Materials

Physical Properties of GP Materials

Wetting angle

Sessile drop method measures the angle of contact between a liquid and solid when a drop of standard saline solution is placed on a contact lens surface in air.

Captive bubble method measures the angle of contact between a gas bubble and a polymer surface when a bubble of air floats up against the underside of a flat polymer surface in standard saline solution. In 1978, the CLMA adopted this method as their standard for determining wetting angles on GP materials.

Wilhelmy plate method wetting angles are not measured directly, but are calculated from force measurements as a function of immersion depth of the lens material in saline solution. Two contact angles, an advancing and a receding angle, for a single lens material can be easily determined.


The difference between these two angles is called the contact angle hysteresis. The receding angle measured by the Wilhelmy plate method has been found to be similar to the contact angle measured by the captive bubble technique. But, both angles are needed to completely describe the wetting properties dynamically.

Boston® Materials

Physical Properties of GP Materials

Wetting angle

Methods of determining contact angles



A word of caution must be raised when attempting to use any of these contact angle methods to predict actual on-eye wetting characteristic of a GP lens. The human tear film contains components (mucin, lipid, lactoferrin, lysozyme, etc.) that significantly contribute to on-eye contact lens wettability.

The drop of conditioning solution applied to a lens is quickly replaced by tear fluid upon insertion of the lens in the eye. Tear components can vary significantly from person to person.

Boston® Materials

Physical Properties of GP Materials

Hardness

Hardness is generally measured in one of two ways, either by the **Rockwell R Hardness method** or the **Shore D Hardness method**.

Both methods measure the relative resistance of a GP material to indentation and provide an empirical hardness value intended primarily for quality control purposes.

No simple correlation exists between hardness determined by either of these two methods and the fundamental properties of the material being tested.

To better understand the resistance of a GP material to mechanical wear and tear, modulus and toughness should be measured

| | PMMA | Boston™ II | BostonES™ | Boston™ IV | Boston™ Equalens™ | Boston EO™ | Boston™ Equalens™ II | Boston XO™ | Boston XO₂™ |
|-------------|------|------------|-----------|------------|----------------------|------------|-------------------------|------------|----------------|
| Rockwell R* | 124 | 119 | 118 | 117 | 117 | 114 | 114 | 112 | 100 |
| Shore D* | 91 | 85 | 85 | 84 | 82 | 83 | 81 | 81 | 78 |

Boston® Materials

Physical Properties of GP Materials

Modulus


(MPa): (Flexural Modulus)

A force (stress) is applied to a lens causing the lens to deform (strain). This deformation is measured until breakage occurs or until the deformation reaches some predetermined point.

The flexural modulus is the ratio of the stress to strain and is a measure of how well a material resists deformation.

This quality relates to the stiffness of the plastic and affects its ability to “mask” astigmatism.

This is also an important factor in determining lens design and thickness.



Boston® Materials

Physical Properties of GP Materials


Toughness

(MNm/m³): (Toughness)

In this test, the lens material sample is flexed until it breaks or reaches some predetermined deformation point. The energy that a lens can absorb before it reaches this point is the toughness of the material.

In GPs, this quality relates to lens handling and durability.

Stress-Strain Testing Fracture Resistance†



| | PMMA | Boston™ II | Boston ES™ | Boston™ IV | Boston™ Equalens™ | Boston EO™ | Boston™ Equalens™ II | Boston XO™ | Boston XO ₂ ™ |
|------------|------|------------|------------|------------|-------------------|------------|----------------------|------------|--------------------------|
| Modulus* | 2432 | 1800 | 1900 | 1600 | 1600 | 1600 | 1300 | 1500 | 1160 |
| Toughness* | 2.5 | 3.0 | 3.4 | 2.8 | 2.8 | 2.6 | 0.8 | 2.6 | 2.7 |

* Data on file

† ISO/Fatt method expressed in cgs units. EOP and Dk/t measurements are approximate. Data on file

Boston® Materials

Physical Properties of GP Materials

Ultraviolet Absorber

Ultraviolet radiation (UVR) is found adjacent to visible light on the electromagnetic spectrum. UVR is emitted in a range of wavelengths and measured in nanometers (one billionth of a meter).


The shortest wavelength contains the most powerful energy and is most harmful.

UVR begins at approximately 100 nm and extends to 400 nm and is split into:

1. **UVC** rays (100–200 nm) are absorbed by earth's atmosphere.
2. **UVB** (280 nm–315 nm) pass through cloud cover that cause tanning and sunburn.
3. **UVA** (316 nm–380 nm) rays are nearest to the visible light spectrum and least dangerous of the three wavelengths.

Use of UVR absorbers in contact lenses reduces the amount of UV radiation that reaches the underlying structures of the eye (cornea, crystalline lens, retina, etc.), but does not offer total protection.

Protective eyewear (sunglasses, goggles, etc.) is still recommended for maximum protection.



Boston® Materials

Physical Properties of GP Materials

Ultraviolet Absorber

The presence of UVR absorbers in GP contact lenses may cause fluorescein pattern detail to be less visible when viewed with the customary white light and cobalt blue filter.

A simple method for enhancing fluorescein pattern viewing is to add a #12 yellow Wratten filter over the front of the slit lamp objectives.

Slit lamp filter kits are available from Bausch + Lomb.



07/11/2019

BAUSCH + LOMB

Boston®
Materials

Boston® Materials


Physical Properties of GP Materials

Plasma Treatment

Plasma Treatment is an exciting advancement in the manufacture of GP lenses made with Boston materials. Lenses are placed in a specialized vacuum chamber and bombarded with oxygen ions through the use of a radio frequency generator. The optimized process effectively removes any remaining residuals from the lens manufacturing process from the surface of the lens without changing the bulk material properties.

This process dramatically reduces the wetting angle of the lenses and may improve their comfort for the wearer.

Percent reduction in contact angle after plasma treatment*



It is important to remember that Boston materials are inherently wettable and do not require plasma treatment to make them wettable.

Therefore, we do not expect Boston plasma treated lenses to become non-wetting after prolonged use and/or time.

Bausch + Lomb, Boston® Materials

The Boston™ lens materials have undergone an evolutionary process that includes the following developments:

- **Increased oxygen** while maintaining good wetting and deposit resistance
- **Increased stability** and durability without compromising corneal physiology
- **Improved lens machining** qualities and yields, without sacrificing clinical performance

A significant advancement in Fluoro Silicone Acrylate (FSA) technology occurred with the introduction of the **AERCOR™ chemical architecture**. This unique polymer chemistry permits to **maintain and increase oxygen delivery while reducing silicon**.

Two of these products are Boston EO™ and Boston ES™.

Boston® Materials

Boston ES, Boston EO

In **Boston ES** and **Boston EO**, the majority of the oxygen permeability is already guaranteed by the oxygen-permeable polymer matrix **Aercor**. As a result, the proportion of silicone and fluorine is much lower than with many other materials.

Not only the oxygen permeability is crucial for the suitability of a lens material for certain fitting cases. Also surface wettability, flexibility, stability, refractive index and specific gravity play an important role.

Boston® Materials Recommendation Overview

| Application/Usage | Boston ES™ AERCOR™ Chemistry | Boston EO™ AERCOR™ Chemistry | Boston XO™ Second Generation FSA* | Boston XO™ Third Generation FSA* |
|--|------------------------------------|------------------------------------|---|--|
| Daily wear | Excellent | Excellent | Excellent | Excellent |
| Planned replacement | Not recommended | Not recommended | Excellent | Excellent |
| Thin designs | Excellent | Excellent | Not recommended | Not recommended |
| Ultrathin designs | Excellent | Not recommended | Not recommended | Not recommended |
| <u>Toric designs</u> | Excellent (Hyperopic toric) | Excellent | Excellent (Hyperopic toric) | Excellent |
| <u>Presbyopic designs</u> | Excellent | Excellent | Excellent | Excellent |
| Corneal rehabilitation (post PMMA wear) | Excellent | Good | Good | Good |
| GPC rehabilitation | Excellent | Excellent | Excellent | Excellent |
| Keratoconus | Excellent | Excellent | Excellent | Excellent |
| <u>Aphakia</u> | Fair | Good | Excellent | Excellent |
| Post corneal surgery fitting | Excellent | Excellent | Excellent | Excellent |
| High oxygen demand corneas | Fair | Excellent | Excellent | Excellent |
| Dry eye | Excellent | Excellent | Excellent | Excellent |

07/11/2019

* Fluoro Silicone Acrylate

BAUSCH + LOMB

Boston®
Materials

Boston® Materials

Boston™ II Material

| | | | |
|------------------------|--|--|--|
| Material Generic Name: | itafocon A | | |
| Material Type: | Silicone acrylate | | |
| Indications: | <ul style="list-style-type: none"> • Myopia • Hyperopia | <ul style="list-style-type: none"> • Astigmatism • Keratoconus | |
| Button Diameter | 12.7mm | | |
| Modality: | Daily wear | | |
| Special Applications:* | <ul style="list-style-type: none"> • Toric (front, back, bi-toric) designs • Identified lipid depositors • Habitual lens breakage | | |
| Special Attributes:* | <ul style="list-style-type: none"> • Fitting characteristics similar to PMMA • Handling tint: blue | | |

| Dk† (ISO/Fatt) | Silicon Content | Refractive Index | Specific Gravity | Rockwell Hardness | Shore D Hardness | Modulus | Toughness | UV Absorber |
|-------------------|--------------------|---------------------|---------------------|----------------------|---------------------|----------|-----------------------|----------------|
| 12 | 10-12% | 1.471 | 1.13 | 119 | 85 | 1800 MPa | 3.0MNm/m ³ | none |

* Data on file

† in cgs units

07/11/2019

BAUSCH + LOMB

Boston®
Materials

Boston ES™ Material

Particularly favorable mechanical properties due to the **Aercor** polymer matrix and the very low silicone content.

Lenses from Boston ES are therefore, as far as technically possible, made **very thin**. In the weak minus range, a **reduction in thickness of 30% is achieved**.

| | | |
|--------------------------------|---|--|
| Material Generic Name: | enfluocon A | |
| Material Type: | AERCOR™ fluoro silicone acrylate | |
| Indications: | <ul style="list-style-type: none"> • Myopia • Hyperopia • Astigmatism • Aphakia | <ul style="list-style-type: none"> • Keratoconus • Post surgical • Irregular corneal conditions • Presbyopia |
| Button Diameter | 12.7mm | |
| Modality: | Daily wear | |
| Special Applications: * | <ul style="list-style-type: none"> • Toric (front, back, bi-toric) designs • Aspheric designs • Presbyopic designs (multifocal/bifocal) • Keratoconus designs | <ul style="list-style-type: none"> • Thin designs (0.12mm @ 3.00 D) • Ultra thin designs (0.10mm @ 3.00 D) • GPC (rehabilitation) • Wetting/lens depositing problems or lens breakage problems |
| Special Attributes: * | <ul style="list-style-type: none"> • Exceptional durability and modulus • Exceptional wetting and deposit resistance | <ul style="list-style-type: none"> • Handling tints: blue, ice blue, green, brown, gray, clear |
| AERCOR™ Chemical Architecture: | <ul style="list-style-type: none"> • Low silicon content • Oxygen-permeable "foundation" | <ul style="list-style-type: none"> • Oxygen-permeable crosslinking |

| Dk † (ISO/Fatt) | Silicon Content | Refractive Index | Specific Gravity | Rockwell Hardness | Shore D Hardness | Modulus | Toughness | UV Absorber |
|-----------------|-----------------|------------------|------------------|-------------------|------------------|----------|------------------------|--------------------------|
| 18 | 5-7% | 1.443 | 1.22 | 118 | 85 | 1900 MPa | 3.4 MNm/m ³ | with/without (blue only) |

* Data on file

† in cgs units

Boston® Materials

Boston™ IV Material

| | | |
|------------------------|--|---|
| Material Generic Name: | itafocon B | |
| Material Type: | Silicone acrylate | |
| Indications: | <ul style="list-style-type: none"> • Myopia • Hyperopia • Astigmatism | <ul style="list-style-type: none"> • Keratoconus • Presbyopia |
| Button Diameter | 12.7mm | |
| Modality: | Daily wear | |
| Special Applications:* | <ul style="list-style-type: none"> • Toric (front, back, bi-toric) • Identified lipid depositors • Rehabilitation of PMMA wearers | |
| Special Attributes:* | <ul style="list-style-type: none"> • Fitting characteristics similar to Boston II • Handling tint: blue | |

| Dk† (ISO/Fatt) | Silicon Content | Refractive Index | Specific Gravity | Rockwell Hardness | Shore D Hardness | Modulus | Toughness | UV Absorber |
|-------------------|--------------------|---------------------|---------------------|----------------------|---------------------|----------|------------------------|----------------|
| 19 | 14-16% | 1.469 | 1.10 | 117 | 84 | 1600 MPa | 2.8 MNm/m ³ | none |

* Data on file

† in cgs units

07/11/2019

BAUSCH + LOMB

Boston®
Materials

Boston® Materials

Boston™ Equalens™ Material

| | | | |
|-------------------------------|--|---|--|
| Material Generic Name: | itafluorococon A | | |
| Material Type: | Fluoro silicone acrylate | | |
| Indications: | <ul style="list-style-type: none"> • Myopia • Hyperopia • Astigmatism | <ul style="list-style-type: none"> • Aphakia • Presbyopia | |
| Button Diameter | 12.7mm | | |
| Modality: | <ul style="list-style-type: none"> • Daily wear • Flexible wear • Extended wear | | |
| Special Applications:* | <ul style="list-style-type: none"> • Variable wearing schedules (FW) • Post-corneal surgery fitting • Handling tint: blue | | |
| Special Attributes:* | <ul style="list-style-type: none"> • Fluorinated polymer for improved oxygen delivery • Contains UV absorber | | |

| Dk† (ISO/Fatt) | Silicon Content | Refractive Index | Specific Gravity | Rockwell Hardness | Shore D Hardness | Modulus | Toughness | UV Absorber |
|-------------------|--------------------|---------------------|---------------------|----------------------|---------------------|----------|------------------------|----------------|
| 47 | 13-15% | 1.439 | 1.19 | 117 | 82 | 1600 MPa | 2.8 MNm/m ³ | with only |

* Data on file

† in cgs units

07/11/2019

BAUSCH + LOMB

Boston®
Materials

Boston® Materials

Boston EO™ Material

| | | |
|--------------------------------|--|---|
| Material Generic Name: | enfluocon B | |
| Material Type: | AERCOR™ fluoro silicone acrylate | |
| Indications: | <ul style="list-style-type: none"> • Myopia • Hyperopia • Astigmatism • Aphakia | <ul style="list-style-type: none"> • Keratoconus • Post surgical • Irregular corneal conditions • Presbyopia |
| Button Diameter | 12.7mm | |
| Modality: | Daily wear | |
| Special Applications:* | <ul style="list-style-type: none"> • Hyperopic toric designs • Keratoconus • Presbyopic designs (multifocal/bifocal) | <ul style="list-style-type: none"> • Thin designs (0.12mm @ 3.00 D) • Corneas with high oxygen demand • GPC rehabilitation • Corneal rehabilitation (after stabilization) |
| Special Attributes:* | <ul style="list-style-type: none"> • Fitting/manufacturing characteristics similar to Boston ES • Excellent wetting/deposit resistance | <ul style="list-style-type: none"> • Handling tints: blue, ice blue, green, brown, gray, ice green, electric blue |
| AERCOR™ Chemical Architecture: | <ul style="list-style-type: none"> • Low silicon content • Oxygen-permeable "foundation" | <ul style="list-style-type: none"> • Oxygen-permeable crosslinking |

| Dk† (ISO/Fatt) | Silicon Content | Refractive Index | Specific Gravity | Rockwell Hardness | Shore D Hardness | Modulus | Toughness | UV Absorber |
|-------------------|--------------------|---------------------|---------------------|----------------------|---------------------|----------|------------------------|----------------|
| 58 | 5-7% | 1.429 | 1.23 | 114 | 83 | 1600 MPa | 2.6 MNm/m ³ | with/without |

* Data on file

† in cgs units

07/11/2019

BAUSCH + LOMB

Boston®
Materials

Boston® Materials

Boston™ Equalens™ II Material

| | | |
|-------------------------------|---|---|
| Material Generic Name: | oprifocon A | |
| Material Type: | Fluoro silicone acrylate | |
| Indications: | <ul style="list-style-type: none"> • Myopia • Hyperopia • Astigmatism | <ul style="list-style-type: none"> • Aphakia • Post surgical • Presbyopia |
| Button Diameters | 12.7mm, 27mm | |
| Modality: | <ul style="list-style-type: none"> • Daily wear • Flexible wear | <ul style="list-style-type: none"> • Extended wear |
| Special Applications:* | <ul style="list-style-type: none"> • Corneas with high oxygen demand • Orthokeratology (ortho-k) • Presbyopic designs (multifocal/bifocal) | <ul style="list-style-type: none"> • Post-corneal surgical fitting • Scleral lens designs |
| Special Attributes:* | <ul style="list-style-type: none"> • High oxygen delivery | <ul style="list-style-type: none"> • Handling tints: blue, clear, green |

| Dk† (ISO/Fatt) | Silicon Content | Refractive Index | Specific Gravity | Rockwell Hardness | Shore D Hardness | Modulus | Toughness | UV Absorber |
|-------------------|--------------------|---------------------|---------------------|----------------------|---------------------|----------|------------------------|----------------|
| 85 | 9-10% | 1.423 | 1.24 | 114 | 81 | 1300 MPa | 0.8 MNm/m ³ | with/without |

* Data on file

† in cgs units

07/11/2019

BAUSCH + LOMB

Boston®
Materials

Boston® Materials

Boston XO™ Material

| | | |
|-------------------------------|--|---|
| Material Generic Name: | hexafocon A | |
| Material Type: | Fluoro silicone acrylate | |
| Indications: | <ul style="list-style-type: none"> • Myopia • Hyperopia • Astigmatism • Aphakia | <ul style="list-style-type: none"> • Keratoconus • Post surgical • Irregular corneal conditions • Presbyopia |
| Button Diameters | 12.7mm, 17mm, 21mm, 26mm | |
| Modality: | Daily wear | |
| Special Applications:* | <ul style="list-style-type: none"> • Post-corneal surgery fitting • Planned replacement programs • Presbyopic designs (multifocal/bifocal) • Corneas with high oxygen demand | <ul style="list-style-type: none"> • Orthokeratology (ortho-k) • Scleral designs • Large diameter lenses for non-diseased eyes |
| Special Attributes:* | <ul style="list-style-type: none"> • High oxygen delivery • Stability equaling that of lower Dk materials | <ul style="list-style-type: none"> • Handling tints: ice blue, violet, blue, green, clear |

| Dk† (ISO/Fatt) | Silicon Content | Refractive Index | Specific Gravity | Rockwell Hardness | Shore D Hardness | Modulus | Toughness | UV Absorber |
|-------------------|--------------------|---------------------|---------------------|----------------------|---------------------|----------|------------------------|----------------|
| 100 | 8-9% | 1.415 | 1.27 | 112 | 81 | 1500 MPa | 2.6 MNm/m ³ | with/without |

* Data on file

† in cgs units

07/11/2019

BAUSCH + LOMB

Boston®
Materials

Boston® Materials

Boston XO2™ Material

Very high DK value of 141 when a high degree of oxygen permeability matters. Then balance with other material properties (refractive index, specific gravity, stability). Especially in plus range, the specific weight of XO² is advantageous. It is 1.19 and therefore belongs to the lighter materials. In case of large diameter lenses and OrthoK, the focus is clearly on oxygen permeability.

| | | |
|------------------------|--|---|
| Material Generic Name: | hexafocon B | |
| Material Type: | Fluoro silicone acrylate | |
| Indications: | <ul style="list-style-type: none"> • Myopia • Hyperopia • Astigmatism • Aphakia | <ul style="list-style-type: none"> • Keratoconus • Post surgical • Irregular corneal conditions • Presbyopia |
| Button Diameters | 12.7mm, 17mm, 21mm, 26mm | |
| Modality: | Daily wear | |
| Special Applications:* | <ul style="list-style-type: none"> • Post-corneal surgery fitting • Planned replacement programs • Presbyopic designs (multifocal/bifocal) • Corneas with high oxygen demand | <ul style="list-style-type: none"> • Orthokeratology (ortho-k) • Scleral designs • Large diameter lenses for non-diseased eyes |
| Special Attributes:* | <ul style="list-style-type: none"> • Hyper-transmissability • Stability equaling that of lower Dk materials | <ul style="list-style-type: none"> • Handling tints: ice blue, violet, blue, green |

| Dk† (ISO/Fatt) | Silicon Content | Refractive Index | Specific Gravity | Rockwell Hardness | Shore D Hardness | Modulus | Toughness | UV Absorber |
|-------------------|--------------------|---------------------|---------------------|----------------------|---------------------|----------|------------------------|----------------|
| 141 | 12-13% | 1.424 | 1.19 | 100 | 78 | 1160 MPa | 2.7 MNm/m ³ | with/without |

* Data on file

† in cgs units

07/11/2019

BAUSCH + LOMB

Boston®
Materials

Boston® Materials

Quantum™ I Material

| | |
|------------------------|--|
| Material Generic Name: | siflufocon A |
| Material Type: | Fluoro silicone acrylate |
| Indications: | <ul style="list-style-type: none">• Myopia• Hyperopia• Astigmatism |
| Button Diameter | 12.7mm |
| Modality: | <ul style="list-style-type: none">• Daily Wear• Flexible Wear |
| Special Applications:* | <ul style="list-style-type: none">• Presbyopic designs (multifocal/bifocal)• Toric lens designs |
| Special Attributes:* | <ul style="list-style-type: none">• Handling tints (I) - ice blue |

| Dk† (ISO/Fatt) | Silicon Content | Refractive Index | Specific Gravity | Rockwell Hardness | Shore D Hardness | Modulus | Toughness (MN m/m ³) | UV Absorber |
|-------------------|--------------------|---------------------|---------------------|----------------------|---------------------|----------|-------------------------------------|----------------|
| 33 | 7-8% | 1.428 | 1.25 | 114 | 84 | 1730 MPa | 2.2 | without |

† in cgs units

Boston® Materials

Quantum™ II Material

| | | |
|-------------------------------|--|---|
| Material Generic Name: | hexafocon A | |
| Material Type: | Fluoro silicone acrylate | |
| Indications: | <ul style="list-style-type: none"> • Myopia • Hyperopia • Astigmatism • Aphakia | |
| Button Diameter | 12.7 mm | |
| Modality: | <ul style="list-style-type: none"> • Daily Wear • Flexible Wear • Extended Wear | |
| Special Applications:* | <ul style="list-style-type: none"> • Post-corneal surgery fitting • Presbyopic designs | <ul style="list-style-type: none"> • Corneas with high oxygen demand • Orthokeratology (ortho-k) |
| Special Attributes:* | <ul style="list-style-type: none"> • Exceptional material stability • Excellent oxygen for all wearing schedules (19% EOP) | <ul style="list-style-type: none"> • Excellent wetting/deposit resistance • Handling tints (1) - ice blue |

| Dk [†] (ISO/Fatt) | Silicon Content | Refractive Index | Specific Gravity | Rockwell Hardness | Shore D Hardness | Modulus | Toughness (MN m/m ³) | UV Absorber |
|-------------------------------|--------------------|---------------------|---------------------|----------------------|---------------------|----------|-------------------------------------|----------------|
| 100 | 8-9% | 1.414 | 1.26 | 112 | 81 | 1415 MPa | 2.1 | without |

† in cgs units

07/11/2019

BAUSCH + LOMB

Boston®
Materials

Boston® Materials

Material Specifications*

| | Boston™ II | Boston ES™ | Boston IV™ | Boston™ Equalens™ | Boston EO™ | Boston™ Equalens™ II | Boston XO™ | Boston XO ₂ ™ | Quantum I | Quantum II |
|---|------------|------------|------------|-------------------|------------|----------------------|------------|--------------------------|-----------|------------|
| Permeability (ISO/Fatt) cgs units † | 12 | 18 | 19 | 47 | 58 | 85 | 100 | 141 | 33 | 100 |
| Rockwell Hardness | 119 | 118 | 117 | 117 | 114 | 114 | 112 | 100 | 114 | 112 |
| Shore D Hardness | 85 | 85 | 84 | 82 | 83 | 81 | 81 | 78 | 84 | 81 |
| Refractive Index | 1.471 | 1.443 | 1.469 | 1.439 | 1.429 | 1.423 | 1.415 | 1.424 | 1.428 | 1.414 |
| Modulus (MPa) | 1800 | 1900 | 1600 | 1600 | 1600 | 1300 | 1500 | 1160 | 1730 | 1420 |
| Toughness (MNm/m ³) | 3.0 | 3.4 | 2.8 | 2.8 | 2.6 | 0.8 | 2.6 | 2.7 | 2.2 | 2.1 |
| Silicon Content | 10-12% | 5-7% | 14-16% | 13-15% | 5-7% | 9-10% | 8-9% | 12-13% | 7-8% | 8-9% |
| Wetting Angle (captive bubble) | 20° | 52° | 17° | 30° | 49° | 30° | 49° | 38° | 48° | 49° |
| Dynamic Contact Angle (advanced/receding) | 58°/57° | 52°/50° | 58°/57° | 59°/56° | 62°/60° | 59°/56° | 59°/58° | 50°/40° | 62°/60° | 67°/66° |
| Specific Gravity | 1.13 | 1.22 | 1.10 | 1.19 | 1.23 | 1.24 | 1.27 | 1.19 | 1.25 | 1.27 |

* All data on file

† $\times 10^{-11} \text{cm}^3 \text{O}_2 (\text{cm}) / [(\text{sec}) (\text{cm}^2) (\text{mm Hg})] @ 35^\circ \text{C}$


07/11/2019

BAUSCH + LOMB

Boston®
Materials

Boston® Materials

Available Standard Blanks - Overview



| Material | Dk (ISO/Fatt) | UV Inhib. | Color | Product No. |
|--|---------------|-----------|---|--|
| Boston II | 12 | no | blue | 1000 B |
| Boston® IV | 19 | no | blue | 1800 B |
| Boston ES® | 18 | yes | blue ice blue green brown gray clear | 2200 B 2200 FI 2200 A 2200 E 2200 G 2200 C |
| | 18 | no | blue | 2211 B |
| Boston® Equalens® | 47 | yes | blue | 1700 B |
| Boston® Equalens® II | 85 | yes | blue | 1500 B |
| | 85 | no | red yellow | 1511 R 1511 Y |
| Boston® Equalens® II, 27 mm | 85 | no | clear | 1531 C |
| Boston EO® | 58 | yes | blue ice blue green ice green brown gray el. Blue | 2600 B 2600 FI 2600 A 2600 AI 2600 E 2600 G 2600 D |
| | 58 | no | blue brown grey | 2611 B 2611 E 2611 G |
| Boston XO® | 100 | yes | blue ice blue green violet | 2300 B 2300 F 2300 A 2300 V |
| | 100 | no | red yellow | 2311 R 2311 Y |
| Boston XO® 17 mm Shouldered | 100 | yes | ice blue | 2317 FS |
| Boston XO® 17 mm pre-cut | 100 | no | clear | 2317 CS |
| Boston XO® 17 mm pre-cut | 100 | yes | ice blue clear | 23F1708003584 23C1708003584 |
| Boston XO® 18,5 mm Shouldered | 100 | yes | ice blue clear | 23F1858253585 23C1858253585 |
| Boston XO® 21 mm Shouldered | 100 | yes | ice blue clear | 2321 FS 2321 CS |
| Boston XO® 26 mm | 100 | yes | clear | 2331 C |
| Boston XO2® | 141 | yes | blue ice blue green violet | 2400 B 2400 F 2400 A 2400 V |
| | 141 | no | blue red yellow | 2411 B 2411 R 2411 Y |
| Boston XO2® 17 mm Shouldered | 141 | yes | ice blue clear | 2417 FS 2417 CS |
| Boston XO2® 18,5 mm | 141 | yes | ice blue clear | 24F1858253585 24C185253585 |
| Boston XO2® 21 mm Shouldered | 141 | yes | ice blue clear | 2421 FS 2421 CS |
| Boston XO2® 25 mm | 141 | yes | clear | 2431 C |
| Boston® Envision base curves | 58 | yes | blue | 7,0-8,3 |
| Boston® Multivision base curves | 18 | yes | blue | 7,0-8,3 |

07/11/2019

BAUSCH + LOMB

Boston®
Materials

Boston® Materials

Boston® Material characteristics in practice

| Boston® Material characteristics | Material | | |
|---|-----------------|-----------|-----------|
| | Boston Equalens | Boston ES | Boston EO |
| Selection aid: | | | |
| Suitability as standard material when no specific needs to be considered | ++ | ++ | + |
| Suitability for higher diopters considering refraction index, Dk and weight (more than +5.00 and -8.00 dpt) | ++ | - | 0 |
| Suitability for extra thin lens designs | 0 | ++ | 0 |
| Suitability if patient has difficult tear film conditions | + | ++ | + |
| Suitability for patients with increased tear lipids and greasy eyelid edges | + | ++ | + |
| For increased protein deposits due to wearing conditions | + | ** | + |
| For daily, very long wearing hours | ++ | + | ++ |
| For occasional overnight wear | + | - | + |
| Suitability for Keratokonus under consideration of special cornea conditions: high irregularity --> stability of material --> Dk | ++ | + | ++ |
| Technical properties: | | | |
| Breakage resitance | + | ++ | + |
| Resistance against signs of wear | + | ++ | + |
| Stability along with corneal astigmatism and rotation symmetric lens design | + | + | + |
| Surface wettability | + | ++ | + |

- ++ very good
- + good
- 0 satisfying
- still sufficient
- not recommended

07/11/2019

BAUSCH + LOMB

Boston®
Materials

Boston® Materials

Enzymatic Cleaning

Boston™ One Step liquid enzymatic cleaner is the easiest way to **remove stubborn protein deposits weekly from GP lenses**. This clear, odorless cleaner works together with the disinfecting step, **right in the lens case**.

The convenience of Boston™ One Step liquid enzymatic cleaner should promote patient **compliance** while increasing lens **wearing comfort** and overall **patient satisfaction** with GP lenses.



- Ingredients:** A sterile, aqueous, solution containing proteolytic enzyme (subtilisin) as the active ingredient, and glycerol. Preservative-free.
- How supplied:** 5ml dropper bottle

| FEATURES | PROVIDES | CLINICAL ATTRIBUTES* |
|--------------------------|--|--|
| Cleaning efficacy | <ul style="list-style-type: none"> Reduces the adherence of and helps remove deposits from lenses | <ul style="list-style-type: none"> Provides clean lenses for clear vision Delivers initial and long term wearing comfort |
| Simple, one step process | <ul style="list-style-type: none"> Convenience, ease of use Promotes patient compliance and satisfaction Eliminates tablets and vials | <ul style="list-style-type: none"> Liquid form eliminates dissolving time Clear odorless formula Clinically proven safe and effective |

Boston® Materials

Slit Lamp Filter Kit

The presence of UV absorbers in some GP materials may cause fluorescein pattern detail to be less visible when viewed with the white light and cobalt blue filter most common in slit lamps.

To aid in the evaluation of fluorescein pattern, these Wratten #12 filters attach unobtrusively to the front of the slit lamp using a small strip of Velcro (included).

With the filter in place, the green glow of the illuminated fluorescein is enhanced for easier viewing.



Boston® Materials

Boston™ Laboratory Lens Cleaner

For laboratory and professional use only
(not approved for patient use).

Boston™ laboratory lens cleaner is a solvent-enhanced formulation that effectively and quickly **removes manufacturing residuals** (pitch, wax, solvent). This cleaner also removes **lipids, body oils, and contaminants (found in personal skin care products) from the surfaces of all GP lenses.**

When used to prepare lenses prior to dispensing, Boston™ laboratory lens cleaner helps prevent or **eliminate lens non-wetting and hazy vision** of GP lenses.



| FEATURES | PROVIDES | CLINICAL ATTRIBUTES |
|-------------------------------------|---|---|
| Combination of powerful surfactants | <ul style="list-style-type: none"> Quick, effective, cleaning of lens surfaces | <ul style="list-style-type: none"> Clean, comfortable GP lenses Helps lens surfaces to wet fast and completely |
| 2-propanol 10% | <ul style="list-style-type: none"> Thorough removal of manufacturing residuals Removes stubborn deposits from difficult-to-clean lenses | <ul style="list-style-type: none"> Ensures excellent initial wetting and good vision Helps to resolve heavy deposit and lens contamination problems |

Boston® Materials

Extended Keratometer Range + Vertex Conversion

Extended Keratometer Range

| Extended Keratometer Range with +1.25 D Lens | | Extended Keratometer Range with -1.00 D Lens | |
|--|----------------|--|----------------|
| Actual Drum Reading | Extended Value | Actual Drum Reading | Extended Value |
| 43.00 D | 5013 D | 36.00 D | 3087 D |
| 43.25 D | 5042 D | 36.25 D | 3109 D |
| 43.50 D | 5072 D | 36.50 D | 3130 D |
| 43.75 D | 5101 D | 36.75 D | 3151 D |
| 44.00 D | 5130 D | 37.00 D | 3173 D |
| 44.25 D | 5159 D | 37.25 D | 3195 D |
| 44.50 D | 5188 D | 37.50 D | 3216 D |
| 44.75 D | 5217 D | 37.75 D | 3237 D |
| 45.00 D | 5246 D | 38.00 D | 3259 D |
| 45.25 D | 5276 D | 38.25 D | 3280 D |
| 45.50 D | 5305 D | 38.50 D | 3302 D |
| 45.75 D | 5334 D | 38.75 D | 3323 D |
| 46.00 D | 5363 D | 39.00 D | 3345 D |
| 46.25 D | 5392 D | 39.25 D | 3366 D |
| 46.50 D | 5421 D | 39.50 D | 3388 D |
| 46.75 D | 5451 D | 39.75 D | 3409 D |
| 47.00 D | 5480 D | 40.00 D | 3430 D |
| 47.25 D | 5509 D | 40.25 D | 3452 D |
| 47.50 D | 5538 D | 40.50 D | 3473 D |
| 47.75 D | 5567 D | 40.75 D | 3495 D |
| 48.00 D | 5596 D | 41.00 D | 3516 D |
| 48.25 D | 5625 D | 41.25 D | 3538 D |
| 48.50 D | 5655 D | 41.50 D | 3559 D |
| 48.75 D | 5684 D | 41.75 D | 3581 D |
| 49.00 D | 5713 D | 42.00 D | 36.02 D |
| 49.25 D | 5742 D | | |
| 49.50 D | 5771 D | | |
| 49.75 D | 58.00 D | | |
| 50.00 D | 58.30 D | | |
| 50.25 D | 58.59 D | | |
| 50.50 D | 58.88 D | | |
| 50.75 D | 59.17 D | | |
| 51.00 D | 59.46 D | | |
| 51.25 D | 59.75 D | | |
| 51.50 D | 60.04 D | | |
| 51.75 D | 60.33 D | | |
| 52.00 D | 60.63 D | | |

Vertex* Conversion

| - | | + |
|--------|-------|--------|
| -3.87 | 4.00 | +4.25 |
| -4.00 | 4.25 | +4.50 |
| -4.25 | 4.50 | +4.75 |
| -4.50 | 4.75 | +5.00 |
| -4.75 | 5.00 | +5.25 |
| -5.00 | 5.25 | +5.62 |
| -5.12 | 5.50 | +5.87 |
| -5.37 | 5.75 | +6.12 |
| -5.62 | 6.00 | +6.50 |
| -5.75 | 6.25 | +6.75 |
| -6.00 | 6.50 | +7.00 |
| -6.25 | 6.75 | +7.37 |
| -6.50 | 7.00 | +7.62 |
| -6.62 | 7.25 | +8.00 |
| -6.87 | 7.50 | +8.25 |
| -7.12 | 7.75 | +8.50 |
| -7.25 | 8.00 | +8.87 |
| -7.50 | 8.25 | +9.12 |
| -7.75 | 8.50 | +9.50 |
| -7.87 | 8.75 | +9.75 |
| -8.12 | 9.00 | +10.12 |
| -8.37 | 9.25 | +10.37 |
| -8.50 | 9.50 | +10.75 |
| -8.75 | 9.75 | +11.00 |
| -8.87 | 10.00 | +11.37 |
| -9.37 | 10.50 | +12.00 |
| -9.75 | 11.00 | +12.75 |
| -10.12 | 11.50 | +13.37 |
| -10.50 | 12.00 | +14.00 |
| -10.87 | 12.50 | +14.75 |
| -11.25 | 13.00 | +15.50 |
| -11.62 | 13.50 | +16.12 |
| -12.00 | 14.00 | +16.75 |
| -12.37 | 14.50 | +17.50 |
| -12.75 | 15.00 | +18.25 |
| -13.00 | 15.50 | +19.00 |
| -13.50 | 16.00 | +19.75 |
| -13.75 | 16.50 | +20.50 |
| -14.12 | 17.00 | +21.50 |
| -14.50 | 17.50 | +22.25 |
| -14.75 | 18.00 | +23.00 |
| -15.12 | 18.50 | +23.75 |
| -15.50 | 19.00 | +24.75 |

07/11/2019

* The distance between the lens and the cornea is 12mm.

BAUSCH + LOMB

Boston®
Materials

Boston® Materials

Dioptr to Radius (mm) Conversion

| Dioptr to Radius Conversion Formula $337.5/D = \text{mm}$ Radius to Dioptr Conversion Formula $337.5/\text{mm} = D$ | | | | | |
|--|---------|---------|--------|---------|--------|
| Dioptr | Radius | Dioptr | Radius | Dioptr | Radius |
| 23.00 D | 14.67mm | 39.00 D | 8.65mm | 49.00 D | 6.89mm |
| 24.00 D | 14.06mm | 39.25 D | 8.60mm | 49.25 D | 6.85mm |
| 25.00 D | 13.50mm | 39.50 D | 8.54mm | 49.50 D | 6.82mm |
| 26.00 D | 12.98mm | 39.75 D | 8.49mm | 49.75 D | 6.78mm |
| 27.00 D | 12.50mm | 40.00 D | 8.44mm | 50.00 D | 6.75mm |
| 28.00 D | 12.05mm | 40.25 D | 8.39mm | 50.25 D | 6.72mm |
| 29.00 D | 11.63mm | 40.50 D | 8.33mm | 50.50 D | 6.68mm |
| 30.00 D | 11.25mm | 40.75 D | 8.28mm | 50.75 D | 6.65mm |
| 31.00 D | 10.88mm | 41.00 D | 8.23mm | 51.00 D | 6.62mm |
| 31.25 D | 10.80mm | 41.25 D | 8.18mm | 51.25 D | 6.58mm |
| 31.50 D | 10.71mm | 41.50 D | 8.13mm | 51.50 D | 6.55mm |
| 31.75 D | 10.63mm | 41.75 D | 8.08mm | 51.75 D | 6.52mm |
| 32.00 D | 10.54mm | 42.00 D | 8.04mm | 52.00 D | 6.49mm |
| 32.25 D | 10.46mm | 42.25 D | 7.99mm | 52.25 D | 6.45mm |
| 32.50 D | 10.38mm | 42.50 D | 7.94mm | 52.50 D | 6.42mm |
| 32.75 D | 10.30mm | 42.75 D | 7.89mm | 52.75 D | 6.39mm |
| 33.00 D | 10.22mm | 43.00 D | 7.85mm | 53.00 D | 6.36mm |
| 33.25 D | 10.15mm | 43.25 D | 7.80mm | 53.25 D | 6.34mm |
| 33.50 D | 10.07mm | 43.50 D | 7.76mm | 53.50 D | 6.31mm |
| 33.75 D | 10.00mm | 43.75 D | 7.71mm | 53.75 D | 6.28mm |
| 34.00 D | 9.92mm | 44.00 D | 7.67mm | 54.00 D | 6.25mm |
| 34.25 D | 9.85mm | 44.25 D | 7.63mm | 54.25 D | 6.22mm |
| 34.50 D | 9.78mm | 44.50 D | 7.58mm | 54.50 D | 6.19mm |
| 34.75 D | 9.71mm | 44.75 D | 7.54mm | 54.75 D | 6.16mm |
| 35.00 D | 9.64mm | 45.00 D | 7.50mm | 55.00 D | 6.13mm |
| 35.25 D | 9.57mm | 45.25 D | 7.46mm | 55.25 D | 6.10mm |
| 35.50 D | 9.50mm | 45.50 D | 7.42mm | 55.50 D | 6.08mm |
| 35.75 D | 9.44mm | 45.75 D | 7.38mm | 55.75 D | 6.05mm |
| 36.00 D | 9.37mm | 46.00 D | 7.34mm | 56.00 D | 6.03mm |
| 36.25 D | 9.31mm | 46.25 D | 7.30mm | 56.25 D | 6.00mm |
| 36.50 D | 9.24mm | 46.50 D | 7.26mm | 56.50 D | 5.97mm |
| 36.75 D | 9.18mm | 46.75 D | 7.22mm | 56.75 D | 5.95mm |
| 37.00 D | 9.12mm | 47.00 D | 7.18mm | 57.00 D | 5.93mm |
| 37.25 D | 9.06mm | 47.25 D | 7.14mm | 57.25 D | 5.90mm |
| 37.50 D | 9.00mm | 47.50 D | 7.11mm | 57.50 D | 5.88mm |
| 37.75 D | 8.94mm | 47.75 D | 7.07mm | 57.75 D | 5.85mm |
| 38.00 D | 8.88mm | 48.00 D | 7.03mm | 58.00 D | 5.83mm |
| 38.25 D | 8.82mm | 48.25 D | 6.99mm | 58.25 D | 5.80mm |
| 38.50 D | 8.76mm | 48.50 D | 6.96mm | 58.50 D | 5.77mm |
| 38.75 D | 8.70mm | 48.75 D | 6.92mm | 58.75 D | 5.75mm |

07/11/2019

BAUSCH + LOMB

Boston®
Materials

Example of customizable Boston® Materials Certificate for successfully fitting specialty lenses

Certificate
No. 17-0001

Boston® Materials -
The Original made by Bausch+Lomb

**Company, Customer
Eye Care Professional Name**

is certified for using original Boston™ Materials –
the right choice for corneal health and exceptional
vision.

- Innovative products for challenging visual conditions
- Specific gas permeable lenses for outstanding oxygen delivery
- Wide range of services

Date Forename Surname

Lab logo(s)

BAUSCH+LOMB
Boston®
Materials

07/11/2019

BAUSCH+LOMB

Boston®
Materials